

YOUR ONE-STOP DIN REPLICATION STATION

SCAN-ONE-PRINT-ONE



Whether you need to add new info to a reprint of an existing Donation Identification Number (DIN) label, or just need to make a copy, this is the perfect solution for you. The Enhanced DIN Replicator is a stand-alone system that includes a barcode printer, scanner and keyboard display. The system is set up to replicate 7-digit Codabar and 13-digit ISBT 128 barcodes by simply scanning an existing DIN. The keyboard allows you to enter in any additional information that isn't available from scanning the old barcode.



- Small, inexpensive, easy to use
- System includes barcode printer, scanner and keyboard display unit
- Replicates Codabar and ISBT 128 DINs by simply scanning an existing label
- Or supplement barcode data with manual input



SAVE LIVES, SAVE TIME, SAVE MONEY

There's a lot of pressure put on the labels you use. If they don't perform well when it comes to readability, trackability and just hanging on, there's more at stake than just a missing label. With adhesives and materials that can withstand harsh environments, and barcode encodation that's been checked and double-checked, Computype's got you covered with products that perform.

LOOK TO US FOR ALL YOUR CRITICAL LABELING SOLUTIONS

We work closely with you to provide the right labels to meet all your requirements. And not just the right labels, but everything for your labeling needs. From labels to printers, and scanners to label replication systems, Computype offers this and everything in between.

WE'D LOVE TO HEAR FROM YOU.

Computype
2285 West County Road C
St. Paul, MN 55113-2567
Phone: 800-328-0852 | 651-633-0633
Fax: 651-633-5580
Email: sales@computype.com
Twitter: @computype

www.computype.com

Computype
Kingswood Business Park, Connaught Road
Kingston Upon Hull
HU7 3AP
United Kingdom
Phone: [44] 1482-835366
Fax: [44] 1482-831152 421-606
Email: contactus@computype.com

TM0310.r4.0615



Computype™