

MOTOROLA DS4800 SERIES SCANNERS

- LIGHT WEIGHT • ERGONOMICALLY FRIENDLY
- SCAN LINEAR (ISBT 128) AND 2D SYMBOLS
- IMAGE CAPTURE DONOR ID CARDS, ETC.



Outfit your staff with a light weight (6 oz.) scanner/imager that can accommodate ISBT 128 symbols, 2D symbols, as well as donor cards and other documents with ease. With standard USB cable, this is truly a plug-and-play input device.

- › The sleek design combines the strength of an arc with a grip and trigger position that maximizes scanning comfort
- › Scan linear or 2D symbols on labels or on a donor's mobile phone
- › Customize the handle with your institution's logo/tag line
- › Handheld or hands-free operation with the optional "Intellistand"
- › Easy and accurate aiming with illumination system that clearly defines the scanning field and offers a small, sharp aiming dot that is visible over the entire working range
- › Capacitive touch trigger with haptic feedback (successful scan can be indicated by vibrations, visual feedback or melodic tones)



Specifications

Physical Dimensions	7.47 in. L. x 2.61 in W. x 1.90 H.
Weight	6 oz. including cable weight of 2 oz.
Power Source	USB VBUS
Color	Alpine White and Midnight Black

Performance

Skew tolerance	86°
Pitch tolerance	82°
Roll tolerance	360°
Motion tolerance	Up to 20 in./second
Light sources	Aiming: 606 nm LED; illumination: 640 nm LED

User environment

Operating temp.	32° F to 104° F [0° C to 40° C]
Storage temp.	-40° F to 158° F [-40° C to 70° C]
Humidity	5% to 95% RH, non-condensing
Drop spec.	Withstands multiple 51 ft. drops to concrete
Tumble testing	250 cycles [500 drops] in ½ meter tumbler

Barcode symbologies

Linear	ISBT 128, incl. concatenated; Code 128, Code 39; Codabar; Code 93; Code 11; others
2D	DataMatrix; PDF417; QR Code; Maxicode, Aztec; others

WE HAVE WHAT YOU NEED FOR SCANNING AND PRINTING

Partner with Computype to enjoy ALL the productivity benefits of barcode scanning: Data collection that's lightning fast and incredibly accurate. We make sure your entire system works together—from the labels to the scanners to the printers, all the way to print-and-apply automation, if that's what's required. We have the tools you need to get the results you want.

LOOK TO US FOR ALL YOUR CRITICAL LABELING SOLUTIONS

We work closely with you to provide the right labels to meet all your requirements. And not just the right labels, but everything for your labeling needs. From labels to printers, and scanners to label replication systems, Computype offers this and everything in between.

WE'D LOVE TO HEAR FROM YOU

Computype
2285 West County Road C
St. Paul, MN 55113-2567
Phone: 800-328-0852 | 651-633-0633
Fax: 651-633-5580
Email: sales@computype.com
Twitter: [@computype](https://twitter.com/computype)

www.computype.com

Computype
Kingswood Business Park, Connaught Road
Kingston Upon Hull
HU7 3AP
United Kingdom
Phone: [44] 1482-835366
Fax: [44] 1482-831152 421-606
Email: contactus@computype.com

www.computype.com

C00007.r0.0914