

# GET YOUR BLANK LABELS ANY WAY YOU NEED THEM.

ISBT 128 BLANK LABEL STOCK

If your label vendor is offering only one material choice, do you really have a choice? Computype, the largest and most trusted label supplier in the blood bank community, has plenty of choices of blank labels in stock and ready to ship. No matter what printer or software you use, we have competitive pricing on multiple materials. Of course, all are compliant with both US FDA 21 CFR 175.105 and relevant ISEGA standards.

Here's a sampling:

- **TS864:** 2.6 mil glossy polypropylene film (-65°F to 200°F), good tear resistance; good initial tack and ultimate adhesion to blood bags and tubes
- **TS933:** 4.2 mil polypropylene film (-65°F to 200°F), great dispensability; good low-temperature performance, resistance to UV light, water and plasticizers
- **TS577:** 3 mil low-cost paper (-65°F to 200°F), good initial tack and adhesion to paper and plastics; short repositioning time to correct misplaced labels



## SAVE LIVES, SAVE TIME, SAVE MONEY

There's a lot of pressure put on the labels you use. If they don't perform well when it comes to readability, trackability and just hanging on, there's more at stake than just a missing label. With adhesives and materials that can withstand harsh environments, and barcode encodation that's been checked and double-checked, Computype's got you covered with products that perform.

## LOOK TO US FOR ALL YOUR CRITICAL LABELING SOLUTIONS

We work closely with you to provide the right labels to meet all your requirements. And not just the right labels, but everything for your labeling needs. From labels to printers, and scanners to label replication systems, Computype offers this and everything in between.

## WE'D LOVE TO HEAR FROM YOU.

Computype  
2285 West County Road C  
St. Paul, MN 55113-2567  
Phone: 800-328-0852 | 651-633-0633  
Fax: 651-633-5580  
Email: [sales@computype.com](mailto:sales@computype.com)  
Twitter: [@computype](https://twitter.com/computype)

[www.computype.com](http://www.computype.com)

Computype  
Kingswood Business Park, Connaught Road  
Kingston Upon Hull  
HU7 3AP  
United Kingdom  
Phone: [44] 1482-835366  
Fax: [44] 1482-831152 421-606  
Email: [contactus@computype.com](mailto:contactus@computype.com)

TM0302.r0.0615