

# Save time and maintain accuracy with pathology specimen labeling.

HISTOLOGY CASSETTE LABELS AND PATHOLOGY SLIDE LABELS

## Histology Cassette Labels

We've engineered a new line of easy-to-use histology cassette labels. Designed for application on cassettes prior to the tissue fixing process, this new label provides adhesion and image stability without necessitating the purchase and maintenance of any additional labeling equipment. These labels are exactly what you would expect from Computype:

- Easy to use — just print-and-apply
- Engineered and tested to withstand tissue fixing process
- Permanently adhered to cassettes with no additional equipment required



## Pathology Slide Labels

Every year, tens of millions of Computype labels get used in pathology labs around the world. Is your lab using them? If not, you're missing out on two decades of experience in creating, providing and improving labeling solutions for pathology labs like yours. It's this understanding that leads to great innovations like our Flexible Lamination After Printing (FLAP) slide label. Uniquely constructed, our pathology FLAP labels are:

- Manufactured to withstand the chemical environments and processes of pathology staining
- Designed specifically for use on microscopic slides and/or samples
- Over-laminated for extra protection, ensuring 100% accuracy and scannability

LS0909.r1.0613

### **FORMIDABLE LAB ENVIRONMENT, MEET FORMING A PARTNERSHIP THAT MAKES THE GRADE.**

From printers to scanners to label replication systems, Computype can equip and provide solutions for you no matter the challenges of your lab environment. Over the past twenty years, we've put our labels and adhesives through endless tests to prove they can keep up with lab demands of your laboratory.

### **WE'D LOVE TO HEAR FROM YOU.**

Computype  
2285 West County Road C  
St. Paul, MN 55113-2567

Twitter: @computype  
Email: sales@computype.com  
Fax: 651-633-5580  
Phone: 800-328-0852 | 651-633-0633  
computype.com