


LabTrack™ Software


LABORATORY SAMPLE MANAGEMENT

GAIN COMPLETE END-TO-END VISIBILITY IN YOUR LAB. LABTRACK™ SOFTWARE EMPOWERS DIGITALLY-CONNECTED LABORATORIES TO EASILY MANAGE SAMPLE INVENTORIES AND A DETAILED CHAIN OF CUSTODY – SAVING TIME, IMPROVING THROUGHPUT AND IMPROVING ACCURACY.



Cloud-based identification and tracking software solution

At its core, LabTrack™ is a comprehensive software solution that enables you to track and monitor individual assets, ensuring a precise chain of custody and real-time visibility into item history, movement, and status. Your LabTrack™ advantage includes:




FLEXIBLE DATA ENTRY

LabTrack™ has three methods for data entry: manual text entry, barcodes, and RFID.



PRINTING & RE-PRINTING

Streamline your workflow by easily printing and re-printing item labels directly from LabTrack™.



RFID & MOBILE CAPABILITIES

With RFID and mobile capabilities, LabTrack™ ensures seamless monitoring and control.




SUSTAINABLE SOLUTIONS

LabTrack™ enables sustainable practices by supporting the recycling and re-use of labware.

Locate, store and monitor samples

LabTrack™ is a dynamic laboratory asset and data management solution designed to streamline operations and elevate efficiency.

- Accurately track item-level sample data
- Instantly pinpoint sample locations
- Conveniently monitor sample status and usage




WE'D LOVE TO HEAR FROM YOU.

Computype
2285 West County Road C
St. Paul, MN 55113-2567


Email: sales@computype.com
Phone: 800-328-0852 | 651-633-0633
computype.com

Simplified LabTrack™ Process Flow

LabTrack™ can help you track the movement of your samples throughout your supply chain and storage as well as your scientific processes. Once an item is associated you will have complete visibility to that item's history and status, automatically generating a chain of custody as it moves. Once your sample reaches the end of its life you can end the journey by deactivating the item ID, or you can recycle the virtual barcode and reassociate it with a new item.

Search 

Search for any specific item at any point in your workflow.



POWERED BY THE ATMA.IO CONNECTED PRODUCT CLOUD FROM AVERY DENNISON – SECURELY BUILT ON MICROSOFT AZURE

LabTrack™ Software is a result of our partnership with Avery Dennison. Our software is powered by their atma.io connected product cloud, built on Microsoft Azure for unparalleled security in enterprise-level cloud computing solutions.

With the power and reliability of Microsoft Azure, you can rest assured that your data is secure and always accessible, knowing that the platform's global network of data centers ensures uninterrupted service with 99.9% guaranteed uptime.

Furthermore, Computype guarantees that your data belongs to you, we never use customer data for any purpose.



atma.io

WE'D LOVE TO HEAR FROM YOU.

Computype
2285 West County Road C
St. Paul, MN 55113-2567

Email: sales@computype.com
Phone: 800-328-0852 | 651-633-0633
computype.com

